

Hand-in-hand work between the Sponsor and the BA lab : a 5-year journey feedback toward the successful use of an IFN- γ ELISpot method supporting the clinical development of a new vaccine candidate

Christine Bain¹, Philippe Moris², Florence Nicolas², Claire Serraz¹, Audrey Pabiou¹, Magali Roche¹


1. Active Biomarkers - a KCAS company | www.active-biomarkers.com | + 33 (0) 4 37 70 87 00 | 60F Avenue Rockefeller – 69008 – Lyon – France
2. Osivax | www.osivax.com | + 33 (0) 9 70 30 13 80 | 70 Rue Saint-Jean-de-Dieu – 69007 – Lyon – France

ELISpot for monitoring cellular responses to a flu vaccine

Context

OVX836 is a recombinant protein-based **universal flu vaccine** - developed by Osivax - targeting the highly conserved influenza **nucleoprotein (NP)**.


Osivax, as the Sponsor, and Active Biomarkers, as the bioanalytical lab, have worked hand-in-hand to develop a method for evaluating the OVX836 **immunogenicity**.



Approach

IFN- γ ELISpot method developed to monitor T-cell responses to OVX836.

Performance of the method adapted to the **Context of Use** of the biomarker throughout clinical development of OVX836.



Generated by the CTL ImmunoSpot® Series 5 Micro Analyzer

1. Establishing pre-requisites for documenting immunogenicity of OVX836 vaccine


A. Ensuring sample integrity

Study needs

High-quality samples for securing sample analysis and reducing variability

Solutions

- Standardization of preparation, cryopreservation & shipping of clinical samples
- SOPs for thawing and ELISpot



B. Fit-for-purpose development & validation of an IFN- γ ELISpot method

Study needs

Method development for assessing NP-specific cellular immune responses – in reference samples

Solutions

- Optimal conditions for: antigen, viral strain, incubation time and peptide concentration
- 4 reference samples with low, medium, high response range

Antigen	Viral strain	Incubation time	Peptide conc.
Recombinant NP	H1N1 <input checked="" type="checkbox"/>	24h <input checked="" type="checkbox"/>	Conc A <input checked="" type="checkbox"/>
Peptide pool <input checked="" type="checkbox"/>	H2N2	48h	Conc B

: retained condition for clinical runs

- Method precision: intra-run CV<30% and inter-run/ inter-operator <50%
- QC suitability samples for monitoring clinical analytical runs


2. Evaluating cellular response in vaccinated subjects

Study needs

Exploratory assessment of cellular NP-specific immune responses in subjects vaccinated with OVX836

Solutions

- PBMC prepared & cryopreserved in serum-free medium according to Sponsor's SOP
- T cell responses evaluated with validated ELISpot method (see 1B)
- Control charts for monitoring suitability of ELISpot method during clinical sample analysis



Withanage K, et al., J Infect Dis. 2022
*SFC/10⁶ cells: spot-forming cells per million

- Result: Expansion of NP-specific T cells, 1 week after each administration of OVX836 followed by a contraction phase 3 weeks later
- T cell responses dramatically impacted by PBMC quality \Rightarrow ACTION : set a viability threshold: flagging samples at 80% viability – setting 70% as the acceptance limit


3. Method performance adjusted to change in biomarker status and long-term use in clinical development

Study needs

Immunogenicity at day 8 as a primary endpoint in phase 2a \Rightarrow Upgrading method validation

Solutions

- Fine-tuning precision, establishing method sensitivity and linearity
- Checking consistency across lots of critical reagents



- Precision of the method confirmed
- Linearity of the method (see above) demonstrated from LLOQ
- Consistency of the method across different lots of critical reagents


4. Exploring vaccine immunogenicity at high doses

Study needs

NP-specific immune responses at day 8 as a primary endpoint to explore high doses of OVX836

Solutions

- Streamlined process for sample preparation and analysis
- Flagging samples with low viability
- Monitoring suitability of validated ELISpot method (see 3)



Leroux-Roels I, Lancet Infect Dis. 2023

Increase in NP-specific T cells by OVX836 at day 8

5. Robustness of the method proven across clinical studies

- Close communication with the Sponsor
- Well-trained & committed analysts
- Sample integrity (viability >70% after overnight resting)
- Strict control of the critical reagents (bridging studies)
- Standardization of the analytical workflow: from sample preparation, through cryopreservation and shipping, to sample analysis by ELISpot

The method proved to be robust throughout the clinical development of OVX836